# steptics®

vision<sub>LP</sub>

## steptics vision LP glass

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## Table of Contents

1. Description.....	11
2. Intended use.....	11
3. Included components.....	11
4. Indications for use and target and patient group.....	12
5. General safety instructions .....	12
6. Warnings  .....	12
7. Product selection.....	13
8. Assembly instructions.....	13
9. Alignment instructions.....	14
10. Use.....	15
11. Maintenance.....	15
12. Report on a serious incident.....	15
13. Disposal .....	15
14. Liability .....	16
15. Conformity .....	16
16. Name and address of the manufacturer.....	16
17. Technical data .....	16

<b>MD</b> Medical device	<b>steptics vision LP</b> Instructions for Use	<b>CE</b> Version 1.1 Status: 29.01.26
<p>Information:</p> <ul style="list-style-type: none"><li>• Please read this document carefully before using the product.</li><li>• Follow the instructions for fitting and prosthesis assembly carefully.</li><li>• Please keep this document in a safe place.</li><li>• Instruct the patient in the safe use of the product.</li><li>• For users with visual impairments, this document is available as a PDF file for enlarged viewing on our website.</li></ul>		


## 1. Description

The steptics vision LP glass is a lightweight and resilient prosthetic foot made from a glass fiber-reinforced composite material. It is designed for users with mobility levels 1 to 3 and offers dynamic energy return combined with high robustness and low weight. The design allows integration into modular prosthetic systems and supports a wide range of everyday movements.

## 2. Intended use

As part of a modular prosthetic system, the prosthetic foot replaces the function of the foot and ankle of the lower extremity. The suitability of the product for a specific prosthesis and for the respective patient is assessed by an orthopedic technician.

## 3. Included components

Nr.	Designation
1	Glass fiber prosthetic foot
2	Foot shell
3	Spectra sock

Figure 1: steptics vision LP glass

## 4. Indications for use and target and patient group

- Loss of lower limb, amputation or deformity
- Body weight between 50 kg and 140 kg (depending on category)
- Mobility level 1-3 (limited to unrestricted outdoor walker)
- No known contraindications

## 5. General safety instructions

- The use of a lower limb prosthesis entails a fundamental risk of falling, which results in an increased risk of injury.
- This product is intended exclusively for the use of a single person. Reuse by other persons is not permitted.
- Assembly and adjustment may only be carried out by a qualified orthopaedic technician.
- Only use the product for its intended purpose.
- Handle the product carefully to avoid mechanical damage.
- Avoid exposing the product to unauthorised ambient conditions.
- Check the product for functionality and damage before and after each use and after contact with unauthorised ambient conditions.
- If the functionality of the product changes, its function is restricted or there are signs of damage or wear that prevent safe use, the product must no longer be used. In such a case, an orthopaedic specialist must be contacted immediately.
- All aspects listed in the instructions for use that contribute to the proper and safe use of the product should be explained to the patient by an orthopaedic technician before use.
- The foot must be used with the supplied foot shell.

## 6. Warnings

- Improper use may result in injury or damage to the product..
- At the end of its service life, the product must no longer be used, as safety and function are no longer guaranteed.
- Do not drill, grind, cut or otherwise modify the product to ensure safe use by the patient

## 7. Product selection

Ensure that the selected foot meets the patient's requirements in terms of weight, foot size, and mobility. The product is recommended for mobility levels 2-3 (limited to unrestricted outdoor walkers). The following table contains a selection matrix of the appropriate category based on weight and foot size:

User weight	50 - 65 kg	66 - 82 kg	83 - 99 kg	100 - 118 kg	119- 140 kg
Size 23					
Size 24					
Size 25					
Size 26					
Size 27					
Size 28					
Size 29					
Category selection	1	2	3	4	5

## 8. Assembly instructions

**Caution:** Sharp edges or splinters can cause cuts.

- The prosthetic foot is compatible with standardised modular connecting elements.
- Assembly should be carried out in the intended anatomical position.
- Follow the torque recommendations for screws as specified by the manufacturer.

### Foot shell & spectra sock

1. Carefully insert the prosthetic foot into the Spectra sock.
2. Remove the cosmetic connection ring from the foot cosmetic.
3. Use the straight end of a shoehorn to gently insert the foot with the Spectra sock into the foot cosmetic.
4. Slowly pull the shoehorn upward to position the foot completely in the foot cosmetic.
5. Replace the cosmetic connection cap on the foot cosmetic.
6. Once correctly aligned, attach the Spectra sock to the prosthesis to protect it from dust and dirt.

## 9. Alignment instructions


Figure 2: Bench alignment

### Bench alignment

- Sagittal: Align the foot so that the build-up line runs centrally through the pyramid adapter (Figure 2: Bench alignment).
- Frontal: Align the foot so that it is parallel to the load line of the prosthetic socket. If necessary, adjust the rotation.
- Connect the prosthetic foot to the prosthetic socket or prosthetic knee joint using the pyramid adapter. Follow the instructions for use of the pyramid adapter.
- For assembly with a prosthetic knee joint: Follow the instructions for the prosthetic knee joint.

### Static alignment

- Make sure that your weight is evenly distributed on both legs.
- Check that the prosthesis is the correct length in relation to the other leg.
- Check the internal/external rotation.

### Dynamic alignment

To ensure optimal gait and smooth rolling motion, the prosthesis must be adapted to the user's requirements. In addition to the alignment of the prosthetic foot, the choice of shoes also influences the walking characteristics of the foot.

- Rolling over the toes does not work/requires additional energy and the forefoot feels too stiff: Increase the dorsiflexion angle of the prosthesis.
- The prosthetic foot does not return enough energy/the forefoot push-off is too low: Increase the plantar flexion angle of the prosthesis.

## 10. Use

### Cleaning and care

Rinse with clear fresh water or a damp cloth and clean with a mild soap. Then dry with a soft cloth.

### Environmental conditions

- The product is waterproof and can therefore be used in a wet or damp environment and immersed in water up to 3 metres deep for a maximum of 1 hour.
- It tolerates contact with fresh water, salt water, chlorinated water and mild soaps.
- Urine and sweat are considered inadmissible environmental conditions and may cause premature wear.
- Clean the product after contact with fresh water, moisture and other liquids, chemicals and after soiling in accordance with the care and cleaning instructions.
- Keep the product away from open flames.
- Temperature range: -10 to +60 °C
- Humidity range: 0 % to 90 %

### Service life

**Prosthetic foot:** This product complies with the requirements of ISO 10328 and has been tested with two million load cycles. Depending on the patient's activity and frequency and intensity of use, this can correspond to a maximum service life of up to three years.

**Foot shell & Spectra sock:** The cosmetics and Spectra sock are wear parts and are subject to normal abrasion. Their durability depends on individual use and the respective conditions of use.

## 11. Maintenance

The product must be checked regularly by a qualified orthopaedic technician for external damage and functional limitations. The maintenance intervals depend on the individual intensity of use. The prosthesis components must undergo a visual and functional inspection after the first 30 days of use.

## 12. Report on a serious incident

Any serious incident involving the product must be reported immediately to the manufacturer and the competent authorities, quoting the serial number. The serial number can be found either on the label on the packaging, underneath the pyramid adapter or in the front area on the underside of the prosthetic foot.

## 13. Disposal

The product and packaging must be disposed of in an environmentally friendly manner in accordance with local or national disposal regulations.

## 14. Liability

The manufacturer accepts liability provided that the product is used in accordance with the information and instructions in this document. The manufacturer is not liable for damage caused by failure to observe these instructions, in particular due to improper use or unauthorised modification of the product.

## 15. Conformity

This product complies with the requirements of Regulation (EU) 2017/745 on medical devices and has been tested in accordance with ISO 10328. The CE Declaration of Conformity can be downloaded from the manufacturer's website.

## 16. Name and address of the manufacturer

steptics GmbH  
Gmunder Str. 35  
81379 Munich, Germany  
[www.steptics.com](http://www.steptics.com)  
[support@steptics.com](mailto:support@steptics.com)

## 17. Technical data

Size [cm]	23	24	25	26	27	28	29
Heel height [mm]	10 ± 5						
Built height [mm] (with foot shell)	62						
Average weight without foot shell [g]	289	300	318	330	398	405	413
Max. user weight [kg]	99	118	140	140	140	140	140
Mobility level	1-3						
Material	Glass fiber reinforced plastic						

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